

SQUARE "NEW ART"

The Champions® Implant Square is a screw-shaped endosseous one-piece implant.. The implants are available in diameters of 3.0/3.5/4.5/5.5 mm and lengths of 6/8/10/12/14/16/18 mm. The implants cover a wide range of indications.

The implant is suitable for the following indications:

Insertion in the long-term healed jaw bone (delayed implantation)


Delayed insertion (6-8 weeks after the tooth extraction)

Immediate implantation (implantation directly after the tooth extraction, under certain conditions)

The Champions® Implant Square can be placed in bone in the maxilla and mandible (D1-D4).

Please note the restrictions of indications in the Instructions for Use.

Material: Titanium grade 4



Diameter	Length	Height	Art. No.	Type
3.0 mm	10 mm	8.0 mm	1015	○
	12 mm	8.0 mm	1040	
	14 mm	8.0 mm	1080	
	16 mm	8.0 mm	1120	
3.5 mm	6 mm	8.1 mm	0900	
	8 mm	8.1 mm	1000	
	10 mm	8.1 mm	1020	
	12 mm	8.1 mm	1050	
	14 mm	8.1 mm	1090	
	16 mm	8.1 mm	1130	
	18 mm	8.1 mm	1160	
4.5 mm	6 mm	8.6 mm	0910	
	8 mm	8.6 mm	1010	
	10 mm	8.6 mm	1030	
	12 mm	8.6 mm	1060	
	14 mm	8.6 mm	1100	
	16 mm	8.6 mm	1140	
5.5 mm	6 mm	9.2 mm	0920	
	8 mm	9.2 mm	1012	
	10 mm	9.2 mm	1035	
	12 mm	9.2 mm	1070	
	14 mm	9.2 mm	1110	
	16 mm	9.2 mm	1150	