


# DRILLING & CONDENSING PROTOCOL FOR CHAMPIONS IMPLANTS




備(鑽)孔順序

champions  implants




1

D1, ø 3.5 mm-Square  
(R)Evolution & "NEW ART"




2

D1, ø 4.0 mm-Square  
(R)Evolution & "NEW ART"




3

D2, ø 3.5 mm-Square  
(R)Evolution & "NEW ART"




4

D2, ø 4.0 mm-Square  
(R)Evolution & "NEW ART"




**5**

D3, Ø 4.5 mm-Square  
(R)Evolution & "NEW ART"




**6**

D4, Ø 5.5 mm-Square  
(R)Evolution & "NEW ART"




**7**

D1 + D2, Ø 2.5 mm  
Ball-Head




**8**

D1 + D2, Ø 3.0 mm  
Ball-Head




**9**

D3 + D4, Ø 4.0 mm  
Ball-Head



**10**

D1 + D2, Ø 3.5  
Square Classic



**11**

D3 + D4, Ø 3.5  
Square Classic

09.05.2017  
Bohr- und Condensierungsprotokoll en  
Version 2017-01



champions  implants

Germany Champions Implants (HongKong) Limited  
Room A1, 7/F, Cheuk Nang Plaza,  
250 Hennessy Road, WanChai, HongKong

info@champions-implants.com.hk

More information about our products:  
[www.champions-implants.com](http://www.champions-implants.com)